

DSA 3500HT unit

Technical data



Voltage: 230 V 50/60 Hz

Power rating: 800 W

Main fuses: 2x8 A

Compressor: 90 l/min. 8 bar

Reservoir capacity: 2 x 2 l

Dimensions HxWxD (cm): 105x24x55

Colour: RAL 9010

Option: mast for surgical lighting

Instruments from left to right (front view):

- * Position 1: LUZZANI syringe
- * Position 2: BIEN AIR MC3 micro motor or cord for turbine
- * Position 3: BIEN AIR MC3 micro motor or cord for turbine
- * Position 4: BIEN AIR MC3 micro motor or LED curing light
- * Position 5: EMS or SATELEC scaler
- * Position 6: EMS or SATELEC scaler